

Distributed by: Southern Surgical Inc.

615-327-2644 • 800-800-9729 • southern surgical inc@gmail.com

PRO*facial*



As easy as One....Two.....Three  
Cleanse....Exfoliate....Infuse

Cold Targeting  
Technology



Hair Removal IPL and Laser

DUALMAX



DIOMAX II



Dermaplus  
ipl



# MicroNeedling


## The Science Behind DermaFrac™

DermaFrac™ offers the best infusion of the best topical for the best results.

Vacuum-assisted occlusion places the skin under tension for a pain-free treatment. More importantly, it assists the needle tips to penetrate uniformly to maximize the delivery of the topical infusion solutions into the papillary dermis.

A precision-controlled hand piece with an insertable cartridge containing a high-grade infusion solution rolls across the skin under vacuum.

The vacuum pulls the skin into contact with hundreds of micro-needles that painlessly penetrate through the epidermis and into the upper dermis (papillary layer).



**MORE INFUSION MEANS MORE RESULTS.**

The micro-channels create pathways for introducing the infusion solution into the area of the skin that is most beneficial to the natural rejuvenating process.

The depth of the needles has been carefully selected so that the channels are created just below the dermal-epidermal (DE) junction and just above capillaries and nerve endings. This ensures no bleeding or pain.

# Radio Frequency



Sculpt, Tighten, Enhance  
Immediate & long-term results  
with no discomfort or downtime

Pollogen **LEGEND+** with VoluDerm™ & TriPollar™  
TECHNOLOGIES

# Skin Analysis



# Fractional Laser Skin Resurfacing



## Emerge™ Fractional Laser

EXPAND, BUILD & COMPLEMENT  
Your Aesthetic Practice

- Easy to use, low cost fractional skin resurfacing laser
- High return on investment and great complement to Botox™ Cosmetic and filler treatments to build an aesthetic practice
- Proven technology from the inventor of microfractional laser resurfacing



# Dermo Electroporation

